

Antepartum Fetal Monitor BT-300

BT-300 is showing you what fetus is saying.

Convenient, portable and easy to use, BT-300 is the best antepartum fetal monitor that meets user's satisfaction with all of these distinctive features:

- Compact and light
- High-Sensitivity from the 1 MHz probe
- Simultaneous external twin FHR with dual pulsed doppler
- Ergonomic designed Event Marker
- Economic with outstanding 2-Year Warranty



Fetal Monitor **BT-300**

Specifications

FHR	<ul style="list-style-type: none"> • Pulsed Doppler • Intensity : 10mW/cm² or less • FHR Accuracy : ± 2% of Range 	<ul style="list-style-type: none"> • Ultrasound Frequency : 1 MHz • FHR Range : 30 ~ 240 bpm
FM	<ul style="list-style-type: none"> • Auto-detection of Twin Fetal Movement 	
UC	<ul style="list-style-type: none"> • External Type • Reference (zero) Control 	<ul style="list-style-type: none"> • Frequency Response : DC ~ 0.5 Hz • Measurement Range : 0 ~ 99 units
Printer	<ul style="list-style-type: none"> • Thermal Array Type • Auto Record Period : Off, 10, 20, 30, 40, 50, 60 min 	<ul style="list-style-type: none"> • Print Speed : 1,2,3 cm/min, high speed • Paper Feeding Function
Display	<ul style="list-style-type: none"> • 7-Segment LED 	<ul style="list-style-type: none"> • 3 Channels(FHR I, FHR II, UC)
Sound	<ul style="list-style-type: none"> • Doppler Sound • Information Sound 	<ul style="list-style-type: none"> • Alarm Sound
Set-up	<ul style="list-style-type: none"> • Alarm Upper / Lower Value • Print Speed, Print Contrast 	<ul style="list-style-type: none"> • Time & Date
Function	<ul style="list-style-type: none"> • Mark Function • Zoom-in Function • FM Graph Print 	<ul style="list-style-type: none"> • FHR II offset Function • Auto Print Function
Power	<ul style="list-style-type: none"> • Input : AC 100 ~ 240 V (50/60 Hz) • Consumption : 80 VA 	<ul style="list-style-type: none"> • Output : DC 18 V, 2.5 A
Dimension	<ul style="list-style-type: none"> • 191 mm(H) × 187 mm(L) × 201 mm(D) 	
Standard Accessory	<ul style="list-style-type: none"> • Ultrasound Doppler Probe 2ea or 1ea(Single) • Event Marker Jack 1ea • Power Adaptor & Cord 1ea • Probe Belt 3ea 	<ul style="list-style-type: none"> • UC Probe 1ea • Print Paper 2ea • Ultrasound Gel 1ea • Operation Manual 1ea
Optional Accessory	<ul style="list-style-type: none"> • AST(Acoustic Stimulator) 1ea • Internal Rechargeable Battery 	<ul style="list-style-type: none"> • Bluetooth Antenna
Warranty	<ul style="list-style-type: none"> • 2 Years (Accessory excluded) 	

PC Interface

BCM-330

(Central Monitoring System)

- Remote control up to 8 monitors by central PC
- Simultaneous Fetal heart sound Play and Record
- Unlimited Multi language support
- Wire(RS-232C) and Wireless(Bluetooth) connection

