

p11

E200X Trauma Stretcher

In order to design completely for this stretcher, in addition to provided flexibility and convenience. To design X-Ray cassette tray specially for whole body, it will let the patient with serious flesh wound to avoid moving and transporting trouble.



Specifications

E200X Trauma Stretcher

- ✦ Integrated storage tray with oxygen bottle trough
- ✦ Trend./Reverse Trend.
- ✦ Control hydraulic device in dual-side foot pedal
- ✦ X Ray translucent
- ✦ Brake pedal control four wheels / steering pedal in two end of bed
- ✦ IV pole inserted hole x5
- ✦ Protective corner bumpers x 4
- ✦ Collapsible siderails
- ✦ Footboard bends the angle $0^{\circ} \sim 40^{\circ}$
- ✦ 4" (10cm) mattress or more
- ✦ 8" diameter wheel
- ✦ Wheel Cover
- ✦ Delivery Push Bar
- ✦ Two stages movable IV pole x 1
- ✦ Two stages permanent folding IV rod x 1
- ✦ Pneumatic assist backrest bends the angle $0^{\circ} \sim 90^{\circ}$
- ✦ Foot auxiliary push bar
- ✦ X Ray cassette tray that can slide through all the stretcher (under the platform)
- ✦ Upright oxygen holder aperture

