

E200 Transport Stretcher



Features

There are two kinds of material available for handgrip of side rail choose : one is made of aluminum, the other is made of PP.

Specifications

Cogent E200 Transport Stretcher

It provide sagile and convenience transmission function, the thick comfortable mattress, collapsible and solid of side-rail, controlled system of hydraulic pump, simple and stable control the height and gradient of platform.

Nimble guide installation of four wheels to provide the operator for labor-saving and convenient transport efficiency.

- + Integrated storage tray with oxygen bottle trough
- + Trend / Reverse Trend
- + Double side mounted foot control hydraulic pump
- + Brake pedal control four wheels / steering pedal in two end of bed
- + IV pole inserted hole x 5
- + 3" (8cm) or more mattress
- + 8" diameter wheel
- + Wheel cover
- + Delivery push bar
- + Two stages movable IV rod x1
- + Two stages permanent folding IV Rod x1
- + Foot auxiliary push bar
- + Collapsible siderails
- + Pneumatic assist backrest bends the angle $0^{\circ} \sim 90^{\circ}$
- + Protective corner bumpers x4
- + Upright oxygen holder aperture

