

p09

E200i Emergency Stretcher



Specifications

This stretcher which is continued the merit of E200.
Add the function of knee crook.
It provide the demand and convenience for patient and staff
when they are doing nursing work.

- ✦ Integrated storage tray with oxygen bottle trough
- ✦ Trend./Reverse Trend.
- ✦ Control oil hydraulic device in both foot pedal
- ✦ Drainage Bag clasp x 2
- ✦ Brake pedal control four wheels / steering pedal in two end of bed
- ✦ IV pole inserted hole x5
- ✦ Protective corner bumpers x 4
- ✦ Collapsible siderails
- ✦ Pneumatic assist backrest
- ✦ 8" diameter wheel
- ✦ Wheel Cover
- ✦ Delivery Push Bar
- ✦ Two stages movable IV pole x 1
- ✦ Two stages permanent folding IV rod x 1
- ✦ Backrest bends the angle $0^{\circ} \sim 90^{\circ}$
- ✦ Footboard bends the angle $0^{\circ} \sim 40^{\circ}$
- ✦ Upright oxygen holder aperture

