

Reliance® Ultrasonic Cleaning Systems

Tabletop Ultrasonic Cleaners



Reliance Tabletop Ultrasonic Cleaners

Reliance Ultrasonic Cleaning Systems - Cleaning Beyond the Surface - In Smaller Sizes

The sophistication of today's surgical instruments requires special care and handling prior to the disinfection and sterilization process – on the surface and beyond. You should not have to sacrifice quality if space is limited.

The Reliance Ultrasonic Tabletop Cleaning Systems are designed for processing smaller or individual surgical instruments, when space is restricted. Additionally, they will help you effectively clean, thoroughly rinse and lubricate delicate surgical instruments, in an easy-to-use tabletop design.

Reliance® Ultrasonic Cleaning Systems

Tabletop Ultrasonic Cleaners

Key Features and Benefits

Cleaning Power - 132 kHz Transducers

The Reliance Ultrasonic Tabletop Cleaners are designed with 132 kHz transducers that provide optimal ultrasonic action and noise reduction. This micromechanical scrubbing action can reach small crevices, irregular surfaces or internal passages, allowing complete removal of film and small particles, while protecting delicate instrument surfaces.

All Stainless Steel Construction, Including Tank and Lid

Easy-to-clean, durable stainless steel construction helps extend life of unit while resisting damage caused by chemicals and cracking from dropped instruments.

Digital Controls

Digital control enables precise control of wash cycle time, temperature, and degas cycle.

Lube Pan (optional)

Lube pans are available for instruments requiring lubrication, which promotes instrument longevity and optimal performance. Lube pans are designed to combine the action of cavitation with the soaking of lubricant, and can be stacked on top of unit while a standard mesh basket is being processed, helping to increase productivity.

Two Sizes Available

- 3 ¼ gallon (12L) – 3 ¼ gallon (12L) – CRT3, 11 ¾ x 9 ½ x 8
 - 5 ¾ gallon (22L) – 5 ¾ gallon (22L) – CRT5, 19 ½ x 11 ½ x 6
- Each model is equipped with a drain hose for easy drainage.

Standard Mesh Basket (included)

Baskets are designed with hinged handles that also serve as basket support on top of the unit, which helps promote safety, ergonomics, and fast drying.



Prolystica® Enzymatic Presoak and Cleaner

Recommended STERIS Instrument Cleaning Chemistries for Optimal Cleaning Results

Description	Quantity	Product Number
Prolystica® Enzymatic Presoak and Cleaner	1 gal/ 4 per case	1C3308
Coverage® Spray TB Plus Disinfectant Cleaner	1 gal/ 4 per case	162908
Hinge-Free® Instrument Lubricant	1 gal/ 4 per case	10310

Specifications: For complete specification, please reference STERIS Tech Data Sheet - SD897

STERIS OFFICES WORLDWIDE

Benelux	32 2 523 2488	Japan	81 3 5777 1580
Canada	800 661 3937	Korea	82 2 517 1517
China	86 21 6137 1166	Latin America	800 884 9550
France	33 0 2 38 70 83 50	Malaysia	6 03 7954 9822
Germany	49 2203 890 6969	Singapore	65 68 41 7677
Greece	30 210 6800848	Spain	34 91 658 5920
India	91 33 2367 5150	Switzerland	41 32 376 0200
Italy	39 22 13 034 1	United Kingdom	44 1256 840400

Document # M3264EN.2008-08, Rev. A
Printed 10/2008, 2,500

©2008 STERIS Corporation.
All rights reserved. Printed in USA.

Technologies to Prevent Infection and Contamination™



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 ■ USA
440-354-2600 ■ 800-548-4873
www.steris.com