

the above
curve

V-PRO[®] 60

Low Temperature Sterilization System

Continuing to transform the standard of low temperature sterilization, the V-PRO 60 Sterilizer delivers the flexibility and efficiency you won't find anywhere else...helping you meet the demands of today's high-performing OR's.



Productivity

Multiple Cycles - Having the right cycle for the right instrument means your medical devices are ready when you need them.



Non Lumen Cycle



12 lb load vs. 10.7 lb load



Flexible Cycle



Process longer scopes,
990mm (single & dual) vs. 850mm (single)!




Lumen Cycle




Process single, dual and
triple channel devices

**Similar footprint but
2X the chamber space.**



V-PRO 60
Sterilizer



Competitor's
Chamber

With twice the usable volume of the competitor's sterilizer, the V-PRO 60 Sterilizer helps to maximize your busy OR's productivity.



Ease-of-Use

The V-PRO 60 Sterilizers' easy to learn, user-friendly design can help simplify your day.




VAPROX® HC Sterilant

59% H₂O₂

Same cup can be used in
V-PRO® maX Sterilizer

Minimal Aborts

Moisture check removes residual moisture
from loads to save time and money



Cycle Prompts

Makes it easy to
match the correct
cycle to your
instruments



USB Download

Gather cycle
data instantly



Flexible Installation

Cart or Counter
Mounted,
depending on your needs






Compatibility


Plasma-free processing is gentle on devices.

Chamber H₂O₂ Concentration (%)



Lower concentration equals less potential to damage devices.

Chamber H₂O₂ Concentration²



V-PRO Device Wizard

Quickly identify the medical devices that have been validated by STERIS for processing in the V-PRO Sterilizers



The V-PRO 60 Sterilizer
by the Numbers:
When Less Equals More.

40% Less

Consumable Cost

24% Less

Service Cost

**30% More
Savings***

Only STERIS provides you with an Unbroken Line of Certainty that protects your hospital from what puts it most at risk. A single, interconnected force of:

- Products
- Support
- Workflow Understanding
- Innovation
- Education

All designed to service you best when they are Connected.
Facing unforeseen challenges together with a single vision.

This is the power of....

A
LINE
OF
1

STERIS SecureCareSM Service Solutions

Choose the service that fits you and your budget. With over 800 technicians worldwide, STERIS technicians are available 24/7 when you need it most.

CS-iQ[®] Sterile Processing Workflow Management Software

includes three modules. As your needs evolve, this module can be upgraded to promote productivity, traceability and compliance across your entire department.



To learn more, contact your STERIS representative, visit www.vpro60.com or call 1.800.548.4873.

¹STERRAD NX Technical Information AD-52972 Rev A

²STERRAD NX Data Sheet AD-52368-01-45

³Savings based on Federal Supply Schedule Price List and current as of 8/6/14.

⁴Over 7 year lifetime, when compared to STERRAD NX.

STERRAD 100 NX and STERRAD NX are registered trademarks of Johnson and Johnson Corporation.

Document # M9914EN. 2015-03, Rev. C

©2015 STERIS Corporation.
All rights reserved. Printed in USA.



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-548-4873
www.steris.com